


Calcium Phosphate Coating on Implants made by Additive Manufacturing Methods





Implants made by additive manufacturing methods display porous surface structures suitable for ingrowth of bone tissue. To promote and accelerate bone ingrowth and osseointegration a bioactive BONIT® (Calcium phosphate) coating can be applied on the structured and porous surface.

BONIT® coating on additive manufacturing made implants

"Non-line-of sight" Technology of BONIT®

Due to the "non line of sight" process technology BONIT® coating is also suitable for porous structures, manufactured by additive manufacturing techniques. Deposition from a solution enables a complete and homogenous coverage of porous surface structures.



Titanium-Plasma-Spray (TPS) Surface on an implant TPS + BONIT® on an implant

Properties

Test criteria BONIT®	Result
Coating thickness	~20 µm
Adhesive strength	≥ 15 MPa
Ca/P ratio	1.1 ± 0.1
Phase composition	≥ 70 % brushite / ≤ 30 % HA

Advantages


- Fine crystalline structure with large free surface
- Complete and even coverage of porous surfaces and complex implant shapes
- Simultaneous bone formation and coating substitution
- Replacement by autologous bone
- Microporosity with high capillarity
- Outstanding biocompatibility
- Thin coating thickness
- Optimal solubility
- Faster and better osseointegration
- "Non-line of sight" process

Why BONIT® coating on your structured implant?

BONIT® coating leads to a faster and better osseointegration of your implant.


More Bone Deposition with BONIT®

Significantly more bone ongrowth after 12 weeks through bioactive BONIT® coating was evaluated in an animal study in minipigs [1]




Proven Gap-bridging Capacity

Accelerated gap bridging by BONIT® coating after 4 and 8 weeks in minipigs (circular gap ~ 1 mm) [2]




High Push-out Values

Extremely High Push out values for BONIT® coated samples [3]




Outstanding capillarity

High capillarity of BONIT® coated samples for uptaking of blood and growth factors



Fast Bone Formation *in vivo*

Bone formation on BONIT® coated hip stem, explanted after six weeks



References

1. Schwarz et al.: The transcortical Göttingen Minipig Model for testing osteointegration of implant surfaces and the influence of a soluble CaP coating on osteointegration. Abstract 5901; 5 World Congress of Biomechanics, Munich 2006 (Citation: Journal of Biomechanics 2006; Vol. 39 Suppl.1, page S10).
2. Schwarz et al.: Effect of a resorbable CaP coating on bone-implant contact and density in a gap model after 4 weeks and 8 weeks. An experimental study in Göttinger Minipigs. 52nd Annual Meeting of the Orthopaedic Research Society. Paper No. 0869.2006.
3. S. Tschöfen: Osseointegration von nanokristallinen Diamantoberflächen, Zirkonoxid und Kalziumphosphat im Schafstibiamodel. Diplomarbeit im Rahmen des Dia-Life-Projektes an der Medizinischen Universität Innsbruck, Universitätsklinik für MKG-Chirurgie

Literature

- Szmukler-Moncler S et al. (2003) From Microroughness to Resorbable Bioactive Coatings. In: Bioimplant interface: improving biomaterials and tissue reactions (Eds.: Ellingsen JE & Lyngstadaas SP), Boca Raton, London, New York, Washington D.C.: CRC Press LLC, pp. 73-100.
- McAfee PC et al. (2003) SB Charité Disc Replacement. Report of 60 Prospective Randomized Cases in a U.S. Cente. J Spinal Disorders & Techniques 16 (4), pp. 424-433.
- Becker et al. (2004) Cellular investigation on electrochemically deposited calcium phosphate composites. J Mater Science, Mater in Medicine 15, pp. 437- 440.
- Schwarz et al. (2006) Effect of a Resorbable CaP Coating on Bone-Implant Contact and Density in a Gap Model after 4 and 8 Weeks. An Experimental Study in Göttinger Minipigs. Paper No. 0869, 52nd Annual Meeting of the Orthopaedic Research Society 2006.
- Reigstad et al. (2007) Improved Bone Ingrowth and Fixation With a Thin Calcium Phosphate Coating Intended for Complete Resorption. J Biomed Mater Res Part B: Appl Biomater 83B, pp. 9-15.
- Reigstad et al. (2011) Different patterns of bone fixation with hydroxyapatite and resorbable CaP coatings in the rabbit tibia at 6, 12, and 52 weeks. J Biomed Mater Res B Appl Biomater. 99 (1), pp. 14-20.
- Bontemps G and Schlüter-Brust K (2012) 10 year survival after unicompartmental knee arthroplasty. Prospective long-term follow-up study. 15th ESSKA. Poster
- Schiegnitz et al. (2012) Vertical Osteoconductive Characteristics of Titanium Implants with Calcium-Phosphate- Coated Surfaces – A Pilot Study in Rabbits. Clinical Implant Dentistry and Related Research 16, pp. 194-201.
- Kämmerer et al. (2016) A biphasic calcium phosphate coating for potential drug delivery affects early osseointegration of titanium implants: Journal of Oral Pathology and Medicine. DOI 10.1111
- Thompson K, Petkov S, Zeiter S, Sprecher CM, Richards RG, Moriarty TF, et al. (2019) Intraoperative loading of calcium phosphate-coated implants with gentamicin prevents experimental Staphylococcus aureus infection in vivo. PLoS ONE 14(2): e0210402. <https://doi.org/10.1371/journal.pone.0210402>

We look forward to talking with you!

DOT GmbH
A Member of Eifeler Holdings
Charles-Darwin-Ring 1a
18059 Rostock

Phone +49 381 40 33 50
Fax: +49 381 40 33 599
info@dot-coating.de
www.dot-coating.de



DOT – coating specialist for orthopedic and dental implants

DOT is one of Europe's leading providers of medical coating solutions for orthopedic and dental implants and instruments and also their cleanroom packaging.

Our comprehensive supply chain concept makes us an ideal medical technology partner. Our activities help restore the health of patients worldwide and thus make a major contribution to the improvement of their quality of life.